

## Important Information Concerning Weight Management

Your scale is a helpful tool for monitoring weight. While not the only measure of weight loss, scales are the most popular method used to gauge weight-loss success. A scale measures the sum of your total body weight, which includes bone, muscle, fat, and water.

### Body weight fluctuates during the day, and from day to day, depending on a variety of factors.

- Salt and carbohydrate intake can affect the body's water retention.
- A large meal adds weight and can cause water retention.
- Dehydration from exercise, illness, or low fluid intake can result in weight loss.
- Muscle is also a factor. Heavy-duty resistance training can build muscle, which can affect weight.
- The menstrual cycle can cause temporary weight gains and losses.

## When You're Losing Weight

It's important not to put too much stock in the exact number on the scale, because it can and will vary. This is especially true when you're trying to lose weight. Small weight gains and losses are common as your body adjusts to fewer kilojoules and more exercise. Your scale is a valuable tool when used to track weight over a period of weeks and months. Be aware that different scales often give different results. The scale at your doctor's office may show one weight, and your scale at home another. So don't get too caught up with a single number.

**While a scale can be a useful tool on your wellness journey, it's not a complete program. Come to a WW Workshop [www.ww.com/au/workshop](http://www.ww.com/au/workshop) or join online to learn how to build healthy eating habits, shift your mindset, become more active and understand how to better manage your weight.**

**In Australia: 13 19 97 or [www.ww.com/au](http://www.ww.com/au)**

**In New Zealand: 0800 009 009 or [www.ww.com/nz](http://www.ww.com/nz)**

## Try Not To Over Do It

We generally recommend weighing yourself no more than once a week as part of an overarching wellness plan. Though it's hard to stick to this pattern, it's the best way to prevent unnecessary concern with normal weight fluctuations, which are unrelated to your weight-loss. Regardless of your weighing habits, it's important that you're consistent.

### Here are 5 guidelines for effectively monitoring your weight – weigh yourself:

- Once a week, on the same day of the week
- At the same time of day
- In the same place
- Wearing the same clothes
- Using the same scale

### To keep track of your weight, record only one number for the week.

**This will give you a clear picture of any trend your weight is following.**

## When You're Maintaining Weight

Fluctuations in day-to-day weight tend to be smaller when weight is stable, so a more frequent weigh-in pattern is helpful. To maintain a stable weight:

- Weigh yourself more than once a week.


- Research shows that people who keep weight within a 2-3kg range are more likely to maintain a healthy body weight over the long term.

## Facts You Should Know

Your WW scale is a precise measuring instrument that is most accurate when weighing a stationary object. To ensure accurate readings, always try to stand on the same area of the scale platform and **DO NOT MOVE**. The minimum weight capacity for this scale is 5kg and the maximum weight capacity is 182kg. Any weights above maximum weight capacity will result in an "Err" message.

Electronic sensors are sensitive. Be careful not to drop or jar the scale. Place it gently onto floor surface, and store where it will be protected from impact. The scale is an electronic instrument and should never be submerged in water. Clean with damp cloth.

When the scale batteries need to be replaced (display shows "LO"), remove the screw from the battery cover on the bottom of the scale and remove the cover. Replace the old batteries with 4 new AAA batteries. Replace the battery cover and tighten the screw again. Remove used batteries promptly.



**Dispose of the old batteries properly. Remove batteries from equipment that is not to be used for an extended period of time. Do not attempt to open the scale or to remove any components. Servicing should be done by qualified technicians only. See warranty information below.**

## Before Using Scale


1. Remove all packaging materials.
2. Remove the isolator tab from the battery compartment – without opening the battery cover at the bottom of the scale, gently pull out the isolator tab. (In case the tab cannot be pulled out, remove the screw from the battery cover and remove the cover. Gently pull out the isolator tab. Ensure the battery is installed correctly. Replace the battery cover and tighten the screw again).
3. Set switch on the base of the scale to measure in pounds (lbs) or kilograms (kg).
4. Set scale on floor. Surface must be level for accurate operation and safety.

## The scale should not be used on carpet.

## Instructions for use

Operate the scale on a secure, even floor.

- Step on the scale and stand still. You will see "000" flashing for a few seconds and then your weight reading will appear.



- Step off the scale and it will shut off in 10 seconds automatically. If you keep standing on the scale, the scale will shut off in 15 seconds automatically.

- 'Err' will display if an error occurs, or if your weight is greater than the maximum capacity of the scale.



- Lo' will display if the batteries need to be replaced.



## IMPORTANT !

- Your WW scale is a precise measuring instrument that is most accurate when weighing a person standing still.
- To ensure accurate readings, always try to stand on the same area of the scale's platform and **DO NOT MOVE**.
- Do not get on and off the scale repeatedly. Even a slight change in your position will cause a different reading.
- Ensure the scale is on a flat, level surface.
- Do not keep moving the scale to a new location, as this too can cause different readings. Try to keep it in the same place.
- It is best to measure yourself using the same scale, as readings can vary between scales due to manufacturing tolerances.
- You will get an 'Err' message an error occurs.
- You will get an 'Err' message if your weight is greater than the maximum capacity of the scale.
- To clean the scale, use a damp cloth without using a cleaning product. Make sure that water does not penetrate inside the appliance.
- Do not place anything on the surface of the scale.
- This scale is designed for home use only and is not suitable for professional purposes.

## Battery Notes:

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Replace all batteries of a set at the same time.
- Do not mix old and new batteries. Do not mix alkaline and standard (carbon-zinc) batteries. Do not use rechargeable (Ni-Cad, Ni-MH, etc.) batteries.
- Clean the battery contacts and also those of the device prior to battery installation.
- Ensure the batteries are installed correctly with regard to polarity (+ and -).
- Remove screw to open battery cover and remove old batteries. Remove batteries from scale that is not to be used for an extended period of time.
- "Lo" will appear on the display if battery is low and needs to be replaced. Remove and replace batteries.

When replacing the batteries, make sure you dispose of the old batteries in accordance with applicable regulations. They should be disposed of in a location specifically designated for that purpose so they can be recycled safely and without posing any risk for the environment. Do not burn or bury them.

When using electrical appliances, especially when children are present, basic safety instructions should always be followed, including the following:

- Polythene bags over product or package may be dangerous. To avoid danger of suffocation, keep this wrapper away from babies and children. This bag is not a toy.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

## SAFETY CAUTIONS

- Do not immerse in water or any other liquid.
- Keep the scale out of the reach of children.
- Do not drop the scale as this may damage the sensors.
- To clean the scale, use a damp cloth without using a cleaning product. Make sure water does not penetrate inside the appliance.
- Do not place anything on the surface of the scale.
- Always inspect the appliance before use for noticeable signs of damage.
- Do not use if damaged, or if the appliance has been dropped.
- In the event of damage, contact the Conair Customer Care Line.
- The battery should be removed for safe disposal at the end of its life.
- When the battery has run down, it should be disposed of at a recycling centre to safeguard the environment. If the battery leaks, avoid touching it and dispose of it at a recycling centre.
- This scale is designed for home use only and is not suitable for professional purposes.

When using electrical appliances, especially when children are present, basic safety instructions should always be followed, including the following:

- Polythene bags over product or package may be dangerous. To avoid danger of suffocation, keep this wrapper away from babies and children. This bag is not a toy.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety. Children should be supervised to ensure that they do not play with the appliance.

## LIMITED TEN YEAR WARRANTY

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. In New Zealand our goods come with a guarantee that cannot be excluded under the Consumer Guarantee Act 1993.

Conair Australia Pty Ltd. will repair or replace (at our option, and excludes battery) your bathroom scale free of charge for one hundred and twenty months from the date of purchase if the appliance is defective in workmanship or materials under normal home use.

This warranty covers working parts that affect the function of the appliance. It does not cover cosmetic deterioration caused by fair wear and tear, damage caused by accident, misuse, shipment or other ordinary household use. Conair shall not be responsible for any incidental, consequential or special damages.

This warranty is available to consumers only. You are a consumer if you are the owner of a WW Scale by Conair that was purchased at a retail outlet for personal, family or household use. This warranty is not available to retailers or other commercial purchasers.

This warranty is an additional benefit and does not affect your legal rights. To obtain service under this warranty, please read the Warranty Leaflet included or contact Conair Customer Service at [aus\\_info@conair.com](mailto:aus_info@conair.com).

**While a scale can be a useful tool on your wellness journey, it's not a complete program. Come to a WW Workshop [www.ww.com/au/workshop](http://www.ww.com/au/workshop) or join online to learn how to build healthy eating habits, shift your mindset, become more active and understand how to better manage your weight.  
In Australia: 13 19 97 or [www.ww.com/au](http://www.ww.com/au)  
In New Zealand: 0800 009 009 or [www.ww.com/nz](http://www.ww.com/nz)**

**WW** is the registered trademark of WW International, Inc. and is used under license by Conair Australia Pty Ltd. ©2019 WW International, Inc. All rights reserved.

**WW ASSURANCE:**

Questions? Comments? Call Conair Customer Service.

Conair Australia Pty. Ltd.

Suite 101, 18 Rodborough Rd

Frenchs Forest NSW 2086 Australia

Australia: 1800 650 263

Conair New Zealand Limited

PO BOX 251159

Pakuranga, Auckland 1706 New Zealand

New Zealand: 0800 266 247

For advice & tips on your new WW scale email us on: [aus\\_info@conair.com](mailto:aus_info@conair.com)



## WW STYLE BODY WEIGHT DIGITAL SCALE

**Congratulations! By purchasing this WW Scale, you've made a positive step in your wellness journey. To get the greatest benefit from your purchase, please read this entire brochure before using your scale.**

**WW170A**

IB-19/519