

OPERATING INSTRUCTIONS

1. A plastic cover has been placed over the plug to protect it during transportation. Remove the plastic cover from the plug before use.
2. Ensure the appliance is in the 'off' position when preparing the foot spa for use.
3. Place the foot spa on the floor and fill with warm water. Never add boiling or steaming water. Make sure you do not fill over the 'MAX' line indicated on the inner wall of the foot spa.
4. Plug the appliance into a suitable power outlet, **ensuring that your feet are not in the foot spa when plugging in.**
5. Before placing your feet in the water, make sure you are seated in a chair. You should never stand in your foot spa.
6. Use the on/off control switch to select the 'bubbles/heat' setting.
The invigorating bubbles are released from the bottom of the foot spa through the aeration jets, while the massage tips on the bottom of the foot spa will stimulate the soles of your feet.

The 'bubbles/heat' setting is designed to help maintain the temperature of the warm water. It does not heat the water to a higher temperature.

7. After use, switch the on/off control to the 'off' position. Once your feet are removed from the water, make sure you dry them first before unplugging the unit.
8. Ensure that the foot spa is unplugged from the power outlet before emptying water. See 'User Maintenance' for further instruction.

Arched Pedicure Massager

To use, simply place the soles of your foot on the massager and apply pressure as desired. Move your foot around for a massaging and cleansing feel.

Splash Guard

This foot spa is equipped with a splash guard at the front of the unit to avoid any water splashing over the side. The splash guard is also useful to rest your feet on when you have finished using the foot spa.

Recommended Massaging Times

Generally, the foot spa should not be used for more than 10-15 minutes at any one time. Treatment may be applied once or twice daily as desired. You should adjust your massage time to suit your needs and personal comfort. If you have any concerns regarding treatment, we recommend consulting your physician.

User Maintenance

1. Ensure that the foot spa is switched off and unplugged from the power outlet before emptying water and cleaning. Empty water from the bottom of the foot spa (opposite side to the on/off switch), into a sink or bath. Never submerge the appliance in water.
2. To clean the foot spa, wipe with a damp cloth. Do not use strong chemical cleaners or solvents, harsh detergents or abrasives. Wipe clean and ensure that the appliance is dry and cool before storing. Store the unit out of reach of children.
3. To store the appliance, loosely wrap the power cord around the base of the unit – placing it into the cord storage. Do not wrap the power cord around the top and bottom of the unit, as this could damage the power cord. Never hang the unit by the power cord. Store away in a box or clean, dry place.

LIMITED ONE YEAR WARRANTY

Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. In New Zealand our goods come with a guarantee that cannot be excluded under the Consumer Guarantee Act 1993.

Conair Australia Pty Ltd. will repair or replace (at our option) your appliance free of charge for twelve months from the date of purchase if the appliance is defective in workmanship or materials. This warranty is an additional benefit and does not affect your legal rights. To obtain service under this warranty, please read the Warranty Leaflet included or contact Conair Customer Service at aus_info@conair.com. To make a warranty claim you must retain your proof of purchase.

©2015
Conair Australia Pty. Ltd.
Suite 101, 18 Rodborough Rd
Frenchs Forest NSW, 2086 Australia
Customer Service Hotline Australia:
1800 650 263

New Zealand: 0800 266 247
Email us at: aus_info@conair.com

Body Benefits is a registered trade mark
of Conair Corporation

CFB5CA

IB-15/370

BUBBLING HYDRO SPA RELAXING FOOT BATH



body benefits®
by CONAIR

INSTRUCTION BOOKLET

For your safety and continued enjoyment of this product,
always read the instruction booklet carefully before using

CFB5CA

IMPORTANT SAFETY INSTRUCTIONS

Please read through these instructions carefully and retain for future reference. This appliance is intended for household use.

When using electrical appliances, especially when children are present, basic safety instructions should always be followed.

Polythene bags over product or package may be dangerous. To avoid danger of suffocation, keep this wrapper away from babies and children. This bag is not a toy.

DANGER: When the unit is used in a bathroom, unplug after use since the proximity of water presents a hazard even when the appliance is switched off.

WARNING:

- Always unplug when not in use, and when filling, emptying, cleaning or moving. All appliances are electrically live even when switched off.
- Do not immerse the appliance in water or other liquids.
- Do not reach for an appliance that has fallen into water. Unplug immediately.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure they do not play with the appliance.



IMPORTANT:

- Use this appliance only for its intended use and as described in this instruction booklet.
- Only use on a hard flat surface. Never place on a soft surface or on any surface where it can fall.
- This appliance is not intended for commercial use.
- Take care not to spill any water onto the control switch.
- Do not add boiling water or steaming water to the foot spa - add warm water that is suited to your preference.
- Do not stand in the unit. Always remain seated during use.
- Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock or injury to persons.
- Do not use any other attachments other than those supplied with this product.
- Never drop the appliance or insert any object into an opening, except feet.
- Never fall asleep or use while sleeping.

- Do not use while bathing or in the shower.
- Do not use outdoors or operate where aerosol (spray) products are being used or where oxygen is being administered.
- Keep the cord away from heated surfaces.
- Do not use the appliance if it is damaged. In the event of damage, discontinue use immediately and contact your dealer for instructions on returning it for examination or exchange. If the supply cord is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified person in order to avoid a hazard.
- Never twist or pull on the power cord.
- After use, loosely wrap the cord around the cord storage at the base of the foot spa. Do not wrap tightly as in time this may cause the cord to fracture.
- Always ensure that the voltage to be used corresponds to the voltage marked on the unit.
- For additional protection, the installation of a residual current device (RCD) with a rated operating current not exceeding 30mA is advisable in the electrical current circuit supplying the bathroom. Ask your installer for advice.

CAUTION:

- This appliance should not be used in the following cases: Tuberculosis, benign and malignant tumors, hemorrhages, inflammations of the skin, phlebitis and thrombosis, open and fresh wounds, bruises, broken skin or varicose veins, undetermined calf pain, or joint pain.
- If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- Do not use on any unexplained pain, or swollen muscles, or following a serious muscular injury, before consulting your doctor.
- Consult your doctor if you have any questions regarding the therapeutic use of this foot spa.

BUBBLING HYDRO SPA RELAXING FOOT BATH

Give your tired feet a genuine moment of well-being!

Enjoy an invigorating massage with water, bubbles and heat to help relieve, comfort and stimulate minor muscular aches and pains.

Convenient splash guard for no mess massaging

Custom spa setting combines bubbles and heat*

Aeration jets release millions of bubbles for an invigorating massage

Arched pedicure massager

Massage tips to stimulate the soles of your feet

Extra-deep design to fully immerse both feet in water



Convenient cord storage on base of foot spa

Non-slip rubber feet

1 year warranty